

UNIQUE FEATURES

- Ergonomic Design
- Easy to use

DESCRIPTION

- Autoclave Portable is Single Walled.
- Made of Stainless Steel 304 Q Heavy Gauge.
- The Gasket is made of Superior quality Silicon rubber which is easily Replaceable.
- Electrically operated supplied with pressure gauge, safety valve cum steam release valve.

TECHNICAL DATA

- Working pressure upto 15 PSI
- Power Supply : 230 Volts A.C. 50 Hz Single Phase

OPTIONAL ACCESSORIES

- Stainless Steel Basket
- Digital Temperature Controller & timer

Chamber	Size
Stainless Steel	12" x 12"
Stainless Steel	12" x 20"



MODEL: LSC-05

AUTOCLAVE VERTICAL

UNIQUE FEATURES

- Easy to use
- Negligible Maintenance

DESCRIPTION

- Double walled construction fully made of Stainless Steel 304 Q.
- Electrically operated Supplied with Two Safety Valves, Pressure gauge, Steam Release Valve, Drain Valve at bottom.
- Fitted with Jointless Silicon Rubber Gasket.
- Foot lifting arrangement to open lid.
- I.S.I approved heating elements.

TECHNICAL DATA

- Working pressure upto 15 PSI
- Power Supply : 230 Volts A.C. 50 Hz Single Phase

OPTIONAL ACCESSORIES

- Digital Temperature controller cum indicator with timer.
- Low water level cut off
- Pressure cut off Switch
- Basket Stainless Steel

Chamber Size	Capacity in Liters (Approx)	Rating K.W.
10" X 18"	25	2
12" X 20"	35	2
14" X 22"	50	3
18" X 24"	95	4
18" X 30"	120	4.5
22" X 30"	178	5
22" X 36"	215	6



MODEL: LSC-10

* Specifications subject to change without notice